



Certificate of Containment Testing

Containment Testing of Rotor FA-45-6-30 (5820 715.103-00) in the Eppendorf Centrifuge 5810R

Report No. 40-10A

Report prepared for: Eppendorf AG, Hamburg, Germany
Issue Date: 19th July 2010

Test Summary

Rotor FA-45-6-30 (5820 715.103-00) was containment tested in the Eppendorf Centrifuge 5810R, using Annex AA of IEC 1010-2-20. The sealed bucket was shown to contain a spill within the centrifuge.

Report Written By

A blue ink signature written over a horizontal dashed line.

Report Authorised By

A blue ink signature written over a horizontal dashed line.