

Measurement of sound power level according to
DIN EN ISO 3744:1995-11 and DIN EN ISO 11204:1996-7

Allround mixer MixMate

Eppendorf AG, 22331 Hamburg, Germany

1. Machine:

- 1.1 Type: Allround mixer MixMate
1.2 Serial Number: 5353 05003
1.3 Year of construction: 2008

2. Specifications:

See documents of the manufacturer

3. Measurement system:

- 3.1 Calibrated sound level meter: Norsonic 118

4. Measurement surface:

- 4.1 Measurement surface: hemisphere with 10 measuring points
4.2 Measurement distance: 1 m

5. Measurement conditions:

- 5.1 Environment: Free field over a reflecting plane
5.2 Environmental correction K_{2A}: 0.9 dB(A)
5.7 accuracy class 2
5.8 Operating conditions: 3000 rpm

6. Measurement surface sound pressure level: 46.1 dB(A)
Sound power level: 54.0 dB(A)

TÜV Nord Umweltschutz GmbH & Co. KG

Große Bahnstraße 31, 22525 Hamburg

Nr. 109SST127 Date of measurement: 25.04.2009



Dipl.-Ing. C. Michalke